

Hemosec®
Titanium Clips



Hemosec®
Clip Applicators

Hemosec®

Titanium Clips and Appliers

Product Description:

The Hemosec Clip Applier is a reusable instrument designed to deliver and close metallic ligating clips. This instrument is used with Hemosec Clips which are produced from titanium. The clips are designed for the ligation of tubular structures where a non absorbable ligating device is indicated.

The ligating clips are positioned around a tubular structure and closed by applying pressure to the handles of the applier.

Distinctive Features:

- ◆ Chevron shape of clips facilitates safe “tips-first” closing action.
- ◆ Serrated inner surface and heart-shaped cross-section of the clips allows secure closing.
- ◆ Diamond shaped grooving for atraumatic secure ligation.
- ◆ Triangular shaped wire for easy squeeze application.
- ◆ Adhesive cartridge backing ensures precise and effortless clip loading.
- ◆ **Control:** Jaws respond directly to the surgeon's hand, allowing complete control over clip closure.
- ◆ **Security:** Titanium clip manufacturing process ensures superior retention and reliability.
- ◆ **Comfort:** Ergonomically designed handle provides comfort and control.
- ◆ Colour-coded clips and corresponding colour-coded appliers for different sizes of clips.

TITANIUM CLIPS

Clip Size	Code	Packaging
Small	HTC S-6	Cartridge of 6 Small Clips; Box of 24 Cartridges
Medium	HTC M-6	Cartridge of 6 Medium Clips; Box of 24 Cartridges
Medium-Large	HTC ML-6	Cartridge of 6 Medium-Large Clips; Box of 18 Cartridges
Large	HTC L-6	Cartridge of 6 Large Clips; Box of 18 Cartridges

CLIP APPLIERS

Codes		Sizes
Open Appliers	Laparoscopic Appliers	for Clip Size
LOA-S	—	Small
LOA-M	LLA-M	Medium
LOA-ML	LLA-ML	Medium-Large
LOA-L	LLA-L	Large

Packaging: 1 Applier per Box