



# Esterlus™

## Braided Coated Polyester

---

### Product Details:

Composition	Polyester
Construction	Braided
Coating	Coated with Silicon to provide for smoother and better handling.
In vitro retention of strength	Esterlus sutures are not subject to degradation or loss of tensile strength.
Tissue Reaction	Esterlus sutures elicit only a mild tissue reaction.
Absorption Time	Non-absorbable
Sterilization	100% Ethylene Oxide
Available in Sizes	5-0, 4-0, 3-0, 2-0, 0, 1, 2, 3, 4, 5
Colour	Green and White
Shelf Life	5 years
Packaging	Box of 6 and 12 sutures

### Distinctive Features:







- ◆ Coated with silicon which helps in easy passage through the tissue with minimal trauma.
- ◆ Excellent strength & handling properties.
- ◆ Dyed for good visibility in the surgical arena.
- ◆ Flexible, elongates to support for optimum knotting.
- ◆ Smooth tie-down properties.






# ESTERLUS™ - Non-Absorbable Sutures

## Braided Coated Polyester

### TAPER POINT

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
13 mm Double Armed 	3/8 Circle	75 cm Green		LS 6890	LS 6891				
13 mm Double Armed 	1/2 Circle	75 cm Green			LS 6892				
		75 cm White			LS 6893				
13 mm Double Armed  Pledget 3mm x 3mm x 1.5mm	1/2 Circle	75 cm Green			LS 6892-PL-3				
		75 cm White			LS 6893-PL-3				
25 mm 	1/2 Circle	100 cm Green				LS 655			
25 mm Double Armed  Pledget 3mm x 3mm x 1.5mm	1/2 Circle	90 cm					LS 6226-PL-3		
26 mm Double Armed 	1/2 Circle	75 cm Green				LS 6832			
		100 cm Green				LS 6552			






### TAPER CUT

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
17.5 mm Double Armed 	1/2 Circle	75 cm 5 Green & 5 White					LS 10B52		
		90 cm Green			LS 6935	LS 6936	LS 6937		
		90 cm White					LS 6917		
17.5 mm Double Armed  Pledget 3mm x 3mm x 1.5mm	1/2 Circle	75 cm 5 Green & 5 White					LS 10B55-PL-3		
		90 cm Green			LS 6936-PL-3	LS 6937-PL-3			
		90 cm White				LS 6917-PL-3			
17.5 mm Double Armed  Pledget 6mm x 3mm x 1.5mm	1/2 Circle	75 cm 5 Green & 5 White					LS 10B55		
		90 cm Green					LS 6937-PL-6		
		90 cm White					LS 6917-PL-6		
		90 cm 4 Green					LS 4B37		



# ESTERLUS™ - Non-Absorbable Sutures

## Braided Coated Polyester

### TAPER CUT

Needle Length	Needle Profile	Suture Length	Suture Size							
			3-0	2-0	0	1	2	3	4	5
26 mm Double Armed 	½ Circle	75 cm 5 Green & 5 White		LS 10B72						
		90 cm Green		LS 6977						
		90 cm White		LS 6987						
26 mm Double Armed  Pledget 3mm x 3mm x 1.5mm	½ Circle	75 cm 5 Green & 5 White		LS 10B77-PL-3						
		90 cm Green		LS 6977 PL-3						
		90 cm White		LS 6987-PL-3						
26 mm Double Armed  Pledget 6mm x 3mm x 1.5mm	½ Circle	75 cm 5 Green & 5 White		LS 10B77						
		90 cm Green		LS 6977-PL-6						
		90 cm White		LS 6987-PL-6						
		90 cm 4 Green		LS 4B77						
45 mm HEAVY 	½ Circle	75 cm 1 Green per foil				LS 4842				
		75 cm 4 Green per foil (multi-strand)				LS 4843				
55 mm HEAVY 	½ Circle	75 cm 1 Green per foil							LS 4845	
		75 cm 4 Green per foil (multi-strand)							LS 4846	

### REVERSE CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size						
			4-0	3-0	2-0	0	1	2	3
23 mm 	½ Circle	75 cm Green				LS 6517	LS 6942		
40 mm 	½ Circle	75 cm Green					LS 6941		