



Mass™

Poly(P-dioxanone) Monofilament


Product Details:

Composition	Poly(P-dioxanone)
Construction	Monofilament
Coating	None
In vitro retention of strength	90% at 4 weeks
Tissue Reaction	Non-antigenic, non-pyrogenic and elicits only a mild tissue reaction during the absorption process.
Absorption Time	Completely absorbed in 180 to 240 days
Sterilization	100% Ethylene Oxide
Available in Sizes	6-0, 5-0, 4-0, 3-0, 2-0, 0, 1
Colour	Violet and Undyed
Shelf Life	5 years
Packaging	Box of 12 sutures

Distinctive Features:











- ◆ Ideal for internal tissues where a long lasting, absorbable suture is preferable.
- ◆ Very high initial tensile strength and excellent pliability.
- ◆ Smooth passage through tissue and lesser reactivity.

Tensile Strength Retention:






MASS™ - Synthetic Absorbable Sutures
Poly(P-dioxanone) Monofilament

TAPER POINT



Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
13 mm 	3/8 Circle	70 cm		LNW 1013					
13 mm Double Armed 	3/8 Circle	70 cm		LNW 1013 DA					
13 mm 	1/2 Circle	45 cm			LNW 9102				
		70 cm		LNW 9513 SA					
13 mm Double Armed 	1/2 Circle	45 cm			LNW 9102 DA				
		70 cm	LNW 9613	LNW 9513					
17 mm 	1/2 Circle	90 cm			LNW 9109		LNW 9108		
20 mm 	1/2 Circle	70 cm			LNW 9304	LNW 9237			
20 mm Double Armed 	1/2 Circle	70 cm			LNW 9304 DA	LNW 9237 DA			
30 mm 	1/2 Circle	70 cm				LNW 9332			
30 mm HEAVY 	1/2 Circle	70 cm						LNW 9254	
31 mm 	1/2 Circle	70 cm					LNW 9133		

MASS™ - Synthetic Absorbable Sutures
Poly(P-dioxanone) Monofilament

TAPER CUT

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
30 mm 	½ Circle	70 cm					LNW 9334		
40 mm 	½ Circle	90 cm							LNW 9218
45 mm 	½ Circle	90 cm							LNW 9354

CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
12 mm 	3/8 Circle Premium	45 cm	LNW 9861						
16 mm 	3/8 Circle Premium	45 cm		LNW 9872	LNW 9873				

