



MISCELLANEOUS PRODUCTS

Miscellaneous Medical Products for Various Surgical Procedures:

LOTUS® UMBILICAL COTTON TAPE

Lotus Sterile Umbilical Cotton Tape is available in pre-cut lengths to tie off the umbilical stump during the catheterization process. It is also used to tie off the umbilical cords of newborn infants. The Umbilical Cotton Tape also finds application in a wide range of procedures, as well as for a variety of veterinary purposes.

Product Details:

Composition	White Woven Cotton
Sterilization	100% Ethylene Oxide
Shelf Life	5 years
Packaging	12 pouches packed in a mono carton

Code	Size	Quantity per Pack	Quantity per Box
LW 276	3 mm x 75 cm Cotton Tape	1 length	12 Foils
LW 276 A	3 mm x 25 cm Cotton Tape	3 lengths	12 Foils

LOTUS® TEMPORARY CARDIAC PACING WIRE

Lotus Temporary Cardiac Pacing Wire is a sterile single use product that provides a conductive connection between an external pacemaker and the myocardium. It is used for treating post-operative cardiac arrhythmia of heart surgery.

Product Details:

Composition	Sheathed Multi-filament Steel Wire with Needles
Sterilization	100% Ethylene Oxide
Shelf Life	5 years
Packaging	12 pouches packed in a mono carton

Code	Needle Description	Wire Type	Quantity per Pack	Quantity per Box
LFEP 13	18 mm ½ Circle Taper Cut and 63 mm Straight Cutting Edge, Breakway Needle	Electrode Blue 60 cm Pacing Wire	1	12 Foils
LFEP 15	26 mm ½ Circle Taper Cut and 63 mm Straight Cutting Edge, Breakway Needle	Electrode Blue 60 cm Pacing Wire	1	12 Foils

LOTUS® BONE WAX

Lotus Bone Wax is a sterile mixture of beeswax and isopropyl palmitate, a wax-softening agent, used to help control bleeding from bone surfaces.

Product Details:

Composition	Bees Wax I.P., Iso Propyl Palmitate U.S.P.
Colour	Opaque
Odour	Waxy Odour
Sterilization	Sterilized by Gamma Radiation
Shelf Life	5 years
Packaging	12 pouches packed in a mono carton

Distinctive Features:

- ◆ Helps local haemostasis by functioning like a mechanical barrier.
- ◆ Bio-chemically inert and minimally resorbable.
- ◆ Compatible with body tissues, does not induce skin irritation or any allergic reactions.

Code	Contents per Pack	Quantity per Pack	Quantity per Box
LW 810	2.5 grams	1	12 Foils

LOTUS® PLEDGETS

Lotus Pledgets are sterile single use products, pre-cut rectangular and square-shaped, for use where the surgeon may consider it advantageous to buffer the suture line.

Product Details:

Composition	PTFE (Polytetrafluoroethylene)
Sterilization	100% Ethylene Oxide
Shelf Life	5 years
Packaging	6 pledgets per pouch; 10 pouches per mono carton

Code	Pledget Size	Quantity per Pouch	Quantity per Box
LCP 10	3 mm x 3 mm x 1.5 mm	6 sterile pledgets per foil	10 foils per box
LCP 10-T	3 mm x 3 mm x 0.5 mm	6 sterile pledgets per foil	10 foils per box
LCP 20	6 mm x 3 mm x 1.5 mm	6 sterile pledgets per foil	10 foils per box