



Monolus™

Poly(Glycolide-co-caprolactone) [Poliglecaprone 25]


Product Details:

Composition	Poly(Glycolide-co-caprolactone) [Poliglecaprone 25]
Construction	Monofilament
Coating	None
In vitro retention of strength	41% at 2 weeks
Tissue Reaction	Non-antigenic, non-pyrogenic and elicits only a mild tissue reaction during the absorption process.
Absorption Time	Completely absorbed in 80 to 120 days
Sterilization	100% Ethylene Oxide
Available in Sizes	6-0, 5-0, 4-0, 3-0, 2-0, 0, 1
Colour	Dyed and Undyed
Shelf Life	5 years
Packaging	Box of 12 sutures

Distinctive Features:





- ◆ High initial tensile strength.
- ◆ Soft and smooth feel that helps in easy passage through the tissue with minimal trauma.
- ◆ Excellent handling properties and high pliability.
- ◆ Ideal suture for use in sub cuticular closure.

Tensile Strength Retention:








MONOLUS™ - Synthetic Absorbable Sutures
Poly(Glycolide-co-caprolactone) [Poliglecaprone 25]



TAPER POINT

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
16 mm 	3/8 Circle	70 cm			LNW 1648				
25 mm 	1/2 Circle	70 cm				LNW 1664	LNW 1665		
30 mm 	1/2 Circle	70 cm					LNW 1641	LNW 1642	
36 mm 	1/2 Circle	70 cm					LNW 1666		

CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
12 mm 	3/8 Circle Premium	45 cm	LNW 833						
16 mm 	3/8 Circle Premium	45 cm		LNW 3204	LNW 3205				
25 mm 	3/8 Circle Premium	70 cm				LNW 1326			
25 mm 	5/8 Circle	70 cm					LNW 1602		
60 mm 	Straight	70 cm				LNW 1650			

REVERSE CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
16 mm 	3/8 Circle	70 cm			LNW 1205				
25 mm 	3/8 Circle	45 cm				LNW 3326 (M)-S			
		70 cm				LNW 3326 (M)			