



Nylus®

Black Monofilament Polyamide

Product Details:







Composition	Nylon
Construction	Monofilament
Coating	None
In vitro retention of strength	Nylus is non-absorbable. However, progressive hydrolysis may cause a gradual loss of up to 20% in its tensile strength each year.
Tissue Reaction	Nylus suture elicits a minimal inflammatory reaction in the tissue, followed by gradual encapsulation of the suture by fibrous connective tissue.
Absorption Time	Non-absorbable
Sterilization	100% Ethylene Oxide
Available in Sizes	11-0, 10-0, 9-0, 8-0, 7-0, 6-0, 5-0, 4-0, 3-0, 2-0, 0, 1
Colour	Black
Shelf Life	5 years
Packaging	Box of 12 sutures

Distinctive Features:




- ◆ Easy passage through tissue because of low frictional co-efficient.
- ◆ Excellent elasticity.
- ◆ Knot pull strength and needle swage strength greater than USP limits.
- ◆ Packed using Tyvek paper for extra security.
- ◆ Superior cutting needles which are excellent for use in plastic surgery.
- ◆ Ideal for use in skin closures.

NYLUS® - Non-Absorbable Sutures Black Monofilament Polyamide

CUTTING




Needle Length	Needle Profile	Suture Length	Suture Size							
			6-0	5-0	4-0	3-0	2-0	0	1	
13 mm 	3/8 Circle Premium	70 cm	LNW 3351							
16 mm 	3/8 Circle	35 cm		LNW 3315						
		70 cm			LNW 3318	LNW 3321				
		90 cm			LNW 3319					
19 mm 	3/8 Circle	45 cm			LNW 3356					
		60 cm			LNW 3357					
		70 cm			LNW 3567					
26 mm 	3/8 Circle	70 cm	LNW 3331							
55 mm 	Straight	100 cm						LNW 3255		
60 mm 	Straight	70 cm					LNW 3388	LNW 3390	LNW 3389	

REVERSE CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size							
			6-0	5-0	4-0	3-0	2-0	0	1	
10 mm 	3/8 Circle	38 cm			LNW 3326					
		70 cm	LNW 3320	LNW 3322						
12 mm 	3/8 Circle	70 cm		LNW 3317						
		90 cm		LNW 3316						
12 mm 	3/8 Circle Premium	70 cm			LNW 3752					

NYLUS® - Non-Absorbable Sutures
Black Monofilament Polyamide

SPATULATED

Needle Length	Needle Profile	Suture Length	Suture Size											
			11-0	10-0	9-0	8-0	7-0	6-0	5-0	4-0	3-0	2-0	0	
6 mm 	3/8 Circle Premium Micro Point	38 cm		LNW 3720										
6 mm Double Armed 	3/8 Circle Micro Point	38 cm			LNW 3714	LNW 3324								
6 mm Double Armed 	3/8 Circle Premium Micro Point	38 cm		LNW 3717										

PRE-CUT SUTURES

Suture Length	Suture Size					
	4-0	3-0	2-0	0	1	2
2 x 35 cm			WN 703	WN 704	WN 705	
2 x 76 cm			WN 903	WN 904	WN 905	WN 906

Packaging : Box of 12