



**Pro-AB™**  
**Partially Absorbable Mesh**

# Pro-AB™

## Partially Absorbable Light Weight Mesh for Hernia Repair

### Product Details:

Composition	Poly (glycolide-co-caprolactone), Polypropylene
Construction	Monofilament
Absorption Time	Partially Absorbable
Sterilization	100% Ethylene Oxide
Shelf Life	5 years
Packaging	1 unit per box

Mesh Size	Code
Ø12 cm (Ø4.7")	LAM 1200
6 x 11 cm (2.3" x 4.3")	LAM 66
7.6 x 15 cm (3" x 6")	LAM 114
10 x 15 cm (4" x 6")	LAM 150
15 x 15 cm (6" x 6")	LAM 225
15 x 20 cm (6" x 8")	LAM 300
20 x 30 cm (8" x 12")	LAM 600
30 x 30 cm (12" x 12")	LAM 900

### Distinctive Features:

- ◆ Light weight mesh (Polypropylene amount < 30 gms/m<sup>2</sup>).
- ◆ Excellent handling and better patient comfort.
- ◆ Good initial stiffness during implantation.
- ◆ Excellent burst strength.
- ◆ Good flexibility after absorption.
- ◆ Better tissue integration due to large pores.
- ◆ Blue visible strips aid orientation and placement.
- ◆ Composed from tried and tested materials over a long period of time.

### Indications:

- ◆ Treatment of inguinal and ventral hernias.