



Prosec®

Skin Staplers

Product Description:

The Prosec Skin Stapler is a sterile, single patient use instrument designed to deliver rectangular, stainless steel staples for routine wound closure. The staples are made of lubricant-coated stainless steel.

The Prosec Skin Stapler places single staples to close surgical incisions. It is supplied sterile and is preloaded for single patient use.

Product Details:

Staple Dimensions	Regular	Wide
Wire Diameter	0.50 mm	0.60 mm
Span	5.9 mm	7.2 mm
Leg Length	3.9 mm	4.9 mm

Packaging

Prosec Skin Stapler	Box of 5
Staple Remover	Box of 20

Distinctive Features:

- ◆ Easy staple extraction.
- ◆ Superior cosmetic results - promotes blood supply and faster scar-free healing.
- ◆ Minimizes staple rotation.
- ◆ Uncompromised visibility, comfort and control.
- ◆ Staples are coated with lubricant.

Indications:

The Prosec Skin Stapler has application for routine skin closure in a wide variety of surgical procedures.

Code	Description
LPSS	Prosec Skin Stapler
LSR	Staple Remover