



Solus 910 Swift™

Braided Coated Short-term Poly(Glycolide-co-L-lactide) [Polyglactin-910]


Product Details:

Composition	Poly(Glycolide-co-L-lactide) [Polyglactin-910]
Construction	Braided
Coating	Polycaprolactone and Calcium Stearate
In vitro retention of strength	61% at 1 week
Tissue Reaction	Non-antigenic, non-pyrogenic and elicits only a mild tissue reaction during the absorption process.
Absorption Time	Completely absorbed in 40 to 60 days
Sterilization	100% Ethylene Oxide
Available in Sizes	6-0, 5-0, 4-0, 3-0, 2-0, 0, 1, 2
Colour	Undyed
Shelf Life	5 years
Packaging	Box of 12 sutures

Distinctive Features:



- ◆ High initial tensile strength.
- ◆ Soft and smooth feel that helps in easy passage through the tissue with minimal trauma.
- ◆ Good knotting properties because of better and modern braiding technology.
- ◆ Excellent handling and smooth tie down.

Tensile Strength Retention:





SOLUS 910 SWIFT™ - Synthetic Absorbable Sutures
Braided Coated Short-term Poly(Glycolide-co-L-lactide) [Polyglactin-910]






TAPER POINT

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
26 mm 	½ Circle	70 cm				LNW 2764			
36 mm 	½ Circle	90 cm						LNW 2036	
		100 cm					LNW 2761 RB		

TAPER CUT







Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
36 mm 	½ Circle	70 cm					LNW 2762 ML		
		90 cm					LNW 2762		
		100 cm					LNW 2761		
		110 cm					LNW 2762-ER	LNW 2763	
40 mm 	½ Circle	70 cm						LNW 2763 ML	
		110 cm						LNW 2763 SP	

CUTTING


Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
12 mm 	¾ Circle Premium	45 cm	LNW 9913	LNW 9915					
16 mm 	¾ Circle	45 cm			LNW 2718	LNW 2719			
16 mm 	¾ Circle Premium	75 cm			LNW 9918	LNW 9919			
21.8 mm 	½ Circle	45 cm				LNW 9935			
22 mm 	½ Circle	70 cm				LNW 2735			

SOLUS 910 SWIFT™ - Synthetic Absorbable Sutures
 Braided Coated Short-term Poly(Glycolide-co-L-lactide) [Polyglactin-910]

REVERSE CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size						
			6-0	5-0	4-0	3-0	2-0	0	1
10 mm 	3/8 Circle	45 cm	LNW 2716	LNW 2717					
26 mm 	3/8 Circle	45 cm				LNW 2732 S			
		70 cm				LNW 2732			
26 mm 	3/8 Circle Premium	75 cm					LNW 9933		
35 mm 	3/8 Circle	60 cm						LNW 2721	
36 mm 	3/8 Circle	60 cm					LNW 2720		
40 mm 	1/2 Circle	60 cm					LNW 2720 SP		

DOUBLE ARMED

Needle Length	Needle Profile	Suture Length	Suture Size							
			6-0	5-0	4-0	3-0	2-0	0	1	
40 mm Double Armed 	1/2 Circle Round Bodied, 1/2 Circle Reverse Cutting	90 cm							LNW 2666	
		140 cm						LNW 2777		
		180 cm						LNW 2777 L		