



# Stelus™

## Monofilament 316 L Stainless Steel

---

### Product Details:

Composition	316 L Stainless Steel
Construction	Monofilament
Coating	None
In vitro retention of strength	Stelus sutures have high tensile strength and exhibit very minimal loss over time.
Tissue Reaction	Stelus sutures show very low tissue reactivity.
Absorption Time	Non-absorbable
Sterilization	100% Ethylene Oxide
Available in Sizes	1, 2, 4, 5, 6, 7
Colour	Metallic Silver
Shelf Life	5 years
Packaging	Box of 12 x 6 = 72 sutures Box of 12 x 4 = 48 sutures Box of 12 x 2 = 24 sutures Box of 12 x 1 = 12 sutures




### Distinctive Features:

- ◆ Stelus sutures are laser welded.




# STELUS™ - Non-Absorbable Sutures Monofilament 316 L Stainless Steel



## BLUNT POINT

Needle Length	Needle Profile	Suture Length	Suture Size							
			0	1	2	3	4	5	6	
40 mm 	½ Circle	45 cm (4 x 45 cm)		LS 661	LS 662					
45 mm 	½ Circle	45 cm (4 x 45 cm)					LS 664	LS 651	LS 649	
		75 cm (2 x 75 cm)						LS 643	LS 645	
45 mm 	½ Circle	75 cm (6 x 75 cm)								LS 756 VS

## TAPER CUT

Needle Length	Needle Profile	Suture Length	Suture Size							
			0	1	2	3	4	5	6	
40 mm 	½ Circle	45 cm (4 x 45 cm)		LS 660	LS 650					

## CONVENTIONAL CUTTING

Needle Length	Needle Profile	Suture Length	Suture Size							
			0	1	2	3	4	5	6	7
48 mm 	½ Circle	45 cm (4 x 45 cm)					LS 652	LS 653	LS 654	LS 654-T
		45 cm (6 x 45 cm)							LS 654-6	
		75 cm (2 x 75 cm)						LS 646	LS 644	
57 mm HEAVY 	¾ Circle	75 cm (2 x 75 cm)						LS 647		