

VIROFIL®
Pneumoperitoneum Filter



ProLap®

Laparoscopic Instruments

Virofil® Single-use Pneumoperitoneum Filter

Product Description:

Virofil® is a single-use smoke, plume and pathogen pneumoperitoneum filter device that fits to the trocar, protecting the surgical team during laparoscopy by filtering the contaminated CO₂ that is released from the abdomen.

Distinctive Features:

- ◆ Tested and certified as per FDA norms.
- ◆ Electrostatically charged 0.22 μ hydrophobic filter membrane blocks 99.99% pathogens.
- ◆ Moisture/fluid trap blocks moisture and plume expelled from the abdomen.
- ◆ Smoke trap ensures clear vision throughout surgery.
- ◆ Universal luer lock fits all trocars in the market.

CODE	DESCRIPTION
LVF22	VIROFIL® Smoke, Plume & Pathogen Pneumoperitoneum Filter with Universal Luer Lock + 3 Luer Caps
LNBT5 VF	NOVA™ 5mm Bladeless Trocar with Optical VIEWFINDER, 100mm + VIROFIL™ Smoke, Plume & Pathogen Pneumoperitoneum Filter with Universal Luer Lock + 3 Luer Caps
LNBT10 VF	NOVA™ 5-10mm Bladeless Trocar with Optical VIEWFINDER, 100mm + VIROFIL™ Smoke, Plume & Pathogen Pneumoperitoneum Filter with Universal Luer Lock + 3 Luer Caps
LNBT12 VF	NOVA™ 5-12mm Bladeless Trocar with Optical VIEWFINDER, 100mm + VIROFIL™ Smoke, Plume & Pathogen Pneumoperitoneum Filter with Universal Luer Lock + 3 Luer Caps

Packaging: Instrument provided sterile in 1 unit per pack, 1 sealed pack per box.

ProLap® Single-use Extra-long Trocars

CODE	DESCRIPTION
LBT10XLOP	Single-use Bladeless Trocar 5-10mm with Optical, 150 mm
LBT12XLOP	Single-use Bladeless Trocar 5-12mm with Optical, 150 mm
LBT15XLOP	Single-use Bladeless Trocar 5-15mm with Optical, 150 mm

Packaging: Instrument provided sterile in 1 unit per pack, 1 sealed pack per box.