


Hemosec®
Titanium Clips



Open Applier

Laparoscopic Applier

Hemosec®
Clip Applicators

Hemosec®

Titanium Clips and Appliers

Product Description:

The Hemosec® Clip Applier is a reusable instrument designed to deliver and close metallic ligating clips. This instrument is used with Hemosec® Clips which are produced from titanium. The clips are designed for the ligation of tubular structures where a non absorbable ligating device is indicated.

The ligating clips are positioned around a tubular structure and closed by applying pressure to the handles of the applier.

Distinctive Features:

- ◆ Chevron shape of clips facilitates safe “tips-first” closing action.
- ◆ Serrated inner surface and heart-shaped cross-section of the clips allows secure closing.
- ◆ Diamond shaped grooving for atraumatic secure ligation.
- ◆ Triangular shaped wire for easy squeeze application.
- ◆ Adhesive cartridge backing ensures precise and effortless clip loading.
- ◆ **Control:** Jaws respond directly to the surgeon's hand, allowing complete control over clip closure.
- ◆ **Security:** Titanium clip manufacturing process ensures superior retention and reliability.
- ◆ **Comfort:** Ergonomically designed handle provides comfort and control.
- ◆ Colour-coded clips and corresponding colour-coded appliers for different sizes of clips.

TITANIUM CLIPS

Clip Size	Code	Packaging
Small	HTC S-6	Cartridge of 6 Small Clips; Box of 24 Cartridges
Medium	HTC M-6	Cartridge of 6 Medium Clips; Box of 24 Cartridges
Medium-Large	HTC ML-6	Cartridge of 6 Medium-Large Clips; Box of 18 Cartridges
Large	HTC L-6	Cartridge of 6 Large Clips; Box of 18 Cartridges

CLIP APPLIERS

Codes		Sizes
Open Appliers	Laparoscopic Appliers	for Clip Size
LOA-S	—	Small
LOA-M	LLA-M	Medium
LOA-ML	LLA-ML	Medium-Large
LOA-L	LLA-L	Large

Packaging: 1 Applier per Box